

# TissueGnostics TissueFAXS Cell Analysis System

TissueGnostics

# **TissueFAXS** The equivalent to flow cytometry in tissue

Microscope-based cell analysis system for cells in tissue sections and smears.

TissueFAXS high-throughput acquisition and automated analysis capabilities permit major time savings.





TissueFAXS<sup>®</sup> is a unique analytical instrument that combines the advantages of state of the art:

- multi-channel microscopy
- automated high resolution imaging with the scientific accuracy of flow cytometry.

Instead of bleaching away in the fridge, data of entire fluorescence specimens is available on mouse click to be discussed with your colleagues and students.

The quantitative and objective analysis tools provided by the TissueQuest<sup>®</sup> software module complete the TissueFAXS performance.TissueFAXS is successfully used in cancer research, developmental biology, pathology, immunology, dermatology, urology, drug development and clinical testing.

- Automated multi-channel image acquisition
- Specimen overview and cellular detail in one interface
- Individual regions of interest (ROI) definable by graphical tools
- Automated single cell detection by patented algorithms
- Objective quantification by control based background subtraction
- Dotplot operations and gating on subpopulations
- Backward and Forward Gating from dotplot to image and image to dotplot
- Statistical data in tables and export in Excel as well as ASCII formats





All events in scattergram							
	x-Mean	y-Mean	Events	No./mm <sup>2</sup>	%		
UL	48.24	56.23	463	2625.16	17.00		
UR	71.73	57.48	103	584.00	4.00		
LL	47.08	21.10	1710	9695.50	61.00		
LR	71.39	24.39	542	3073.08	19.00		
SUM	52.85	28.83	2818	15977.73	100.00		

www.tissuegnostics.com

TissueGnostics

# **TissueFAXS** Workflow

Rapid preview	Automated ROI or whole slide acquisition at any magnification	Tunable auto- focus for up to 100x oil objectives	Up to 10 different fluorescence filter cubes	Automatic stitching function generates an exportable sample overview	ROI creation and export

#### Re-acquisition of single out of focus images

#### TissueQuest analysis



www.tissuegnostics.com

PAGE 5

# **TissueFAXS** Management and Acquisition Software

The TissueFAXS acquisition and management software combines an intuitive and user friendly interface and a comfortable workflow with extensive control functions:



- Objective selection
- Light source control
- Filter selection
- Graphical navigation function for 8 slides
- Preview window
- **ROI** management
  - Well plate layout

#### Four colour staining overlay



Nuclei

**Epithelial cells** 

**Basal EC** 



- Live image window
- Stage controls and position info
- Autofocus
- White balance
- Image export
- Colour overlay controls
- TMA layout



Tumor marker

Overlay



#### **TissueGnostics**

# **TissueQuest** Different mask options and quantitative analysis

Analysis of staining intensities can be restricted to the area of the nucleus or can be extended into the cytoplasm. In case of cytoplasmic antigens TissueQuest provides an image processing algorithm that is growing in stained areas around the nucleus.

#### **Ring mask**



Growing mask



# Growing mask/ring mask without nucleus



larkeı

**TissueQuest** Cell Analysis Software Module

**TissueQuest is the first** software solution worldwide to identify and analyse individual cells in their tissue environment

Cell identification is not confined to isolated cells. The TissueQuest software also identifies cells that are agglomerated in dense clusters or have been reduced to ring structures (as it is typical of cancer cells).

Regardless of the type and density of cells involved, TissueQuest recognizes their nuclei due to patented image processing algorithms.

By combining specific stainings and algorithms, one can measure staining intensities that are specific for various cell types.

Cells are clearly distinguished even

- in solid tissue structures
- if embedded in clusters
- if nuclei are uniformly stained
- if staining involves ring structures





48%<sup>+</sup> cells Mean intensity: 47



42%<sup>+</sup> cells Mean intensity: 47



45%<sup>+</sup> cells Mean intensity: 51

#### **TissueGnostics**

# **TissueQuest** Leucocytic infiltrate in rejected kidney: Triple staining

Tissue: Paraffin section of rejected kidney allograft Staining protocol: DAPI, CD68, CD3

DAPI





Cut-offs and/or gates will immediately update the data in the tables.

The export function also outputs the quadrant statistics for all dotplots into Excel and ASCII format.

### **TissueQuest** Quantitative analysis of one marker

Effect of therapy on the density of antigen expression in human skin:

antibody

Control

DAPI

#### Normal control



Normal skin before therapy



Normal skin after therapy



#### By courtesy of Univ. Prof. Dr. Lajos KEMENY, Department of Dermatology and Allergology, University of Szeged, Hungary

It was demonstrated that therapy goes hand in hand with an up-regulation of antigen expression indicated by an increase in the percentage of positive cells from 39% to 45% and in the mean intensity from 39 to 51, respectively.

#### Marker expression before therapy

Negative Control Antibody

MAN AND AND AND

0.5%

99%



Before therapy the percentage of cells that react with the marker was 34% and the mean intensity was 39.

Mean intensity: 29

#### Marker expression after therapy



After therapy the percentage of cells that react with the marker was 45 and the mean intensity was 51.



By courtesy of Ao. Univ. Prof. Dr. Heinz REGELE, Department of Pathology, Medical University of Vienna, Austria

		All events li	scalleryra		
	x-Mean	y-Mean	Events	No./mm <sup>2</sup>	%
UL	48.24	56.23	463	2625.16	17.00
UR	71.73	57.48	103	584.00	4.00
LL	47.08	21.10	1710	9695.50	61.00
LR	71.39	24.39	542	3073.08	19.00
SUM	52.85	28.83	2818	15977.73	100.00

#### www.tissuegnostics.com

#### **TissueGnostics**

#### Tissue staining protocol: DAPI, Plakoglobin-Alexa488, Ki67-APC or IgG<sub>1</sub>-APC

DAPI





Tissue image overlay

Gate1

Gating on

Plakoglobin

DAP

Plakoglobin+cells

60% By the use of gates

defined phenotype.



Cell subpopulations can be visualised in distinctive colours by gating.

# **TissueQuest** Skin Project: Triple staining

The staining combination includes anti-Plakoglobin to identify all keratinocytes and Ki67 to analyse the proliferative activity of all cells. That means that Plakoglobin<sup>+</sup>/Ki67<sup>+</sup> cells are proliferative keratinocytes and Plakoglobin<sup>-</sup>/Ki67<sup>+</sup> cells are proliferative non-epithelial cells.



The use of isotype matched non-reactive control antibodies allows the setting of precise cut-offs.



Statistical results are displayed when cut-off or gates are set. The parameters include mean intensity, relative/absolute events, and percentage of reactive cells as well as sums.

#### 7.15% Plakoglobin+/Ki67+cells

All events in scattergram						
	x-Mean	y-Mean	Events	No./mm <sup>2</sup>	%	
UL	36.81	62.84	192	1451	7.38	
UR	73.78	60.17	186	1406	7.15 ┥	
LL	36.35	0.43	827	6251	31.80	
LR	94.65	0.38	1396	10553	53.67	
SUM	70.35	9.28	2601	19663	100.00	

The quadrant statistics (UL-Upper Left, UR-Upper Right, LL-Lower Left, LR-Lower Right) are presented in tabular format.

#### Plakoglobin

Ki67-APC



By courtesy of Ao. Univ. Prof. Gero KRAMER, Department of Urology, Medical University of Vienna, Austria

By the use of gates the analysis is restricted to subpopulations with



Gates can be used as input gates for new scattergrams

#### TissueGnostics

Tissue: Paraffin section of carcinoma of the prostate Staining protocol: DAPI, basal epithelial cells (basal EC), cytokeratin-18 (prostatic EC), AMACR (tumor cells)







## **TissueQuest** Prostate Project: Four colour staining



Prostatic EC





Analysis of coexpression pattern reveals that 17% of all cells express AMACR and that 99% of the AMACR<sup>+</sup> tumor cells co-express the marker of prostatic EC. In contrast, only 1% of AMACR<sup>+</sup> tumor cells were shown to co-express the marker for basal EC.



AMACR Cytokeratin-18

Statistic for Gate1							
	x-Mean	y-Mean	Events	No./mm <sup>2</sup>	%		
IL	0.32	21.81	568	3680.56	6.00		
IR	86.60	40.60	1707	11061.13	17.00		
L	0.66	9.10	4525	29321.40	48.00		
.R	74.56	10.25	2695	17463.24	29.00		
SUM	37.32	15.85	9495	61526.33	100.00		

#### www.tissuegnostics.com

**PAGE 15** 

Dotplot

# **TissueQuest**

Backward and Forward Gating of single cells & gates

#### **TissueQuest offers an innovative function** for visual and measurement control.

Backward Gating enables the user to connect backwards from a dot in the dotplot to the specific corresponding cell in the image.

By Forward Gating the user is able to locate the dot in the plot corresponding to a specific cell in the image.

The function "View event data" reads out all measurement data for an individual cell.



Overlay



Backward Gating is an easy to use tool that helps in setting cut-offs and in the re-confirmation of image processing data.

It also helps in tracking phenotypically and/or morphologically distinct subsets of cells from the dotplot to the tissue.



# From dot to cell

Image

#### Gate 1



Cytokeratin-18<sup>+</sup> and basal marker negative secretory epithelial cells

#### Gate 2



Cytokeratin-18<sup>+</sup> and basal marker positive basal epithelial cells

#### Gate 3



Cytokeratin-18 negative and basal marker negative non-epithelial cells

www.tissuegnostics.com PAGE 17

Cy3

DAPI

# **TissueFAXS** TMA Module

The TissueFAXS TMA module :

- is located in the new project wizard
- either recognizes cores in full auto mode or provides software-assisted manual core definition
- is tolerant of missing cores
- creates a meta annotation including all properties and a thumbnail of the core
- provides an option for generation of core ID
- allows for the re-naming and change of core location and grouping
- imports and exports from and into Excel sheet templates
- permits automatic acquisition at chosen magnification
- allows for logical grouping
- exports images into TissueQuest and HistoQuest
- provides a TMA-Explorer
- offers a report generator

# 

#### TMA layout

			10 T
	00000	00000	
00000	00000	00000	
06688	00000		
	00000	00000	
00000			
89888		00000	
90060		20000	
88858	00000		
00869		64460	
00000		00000	
00000		00000	
00000		00000	
00000		00000	
09990	1. A.	and the second second	

By courtesy of Univ. Prof. Dr. Kurt ZATLOUKAL, Department of Pathology, Medical University of Graz, Austria







#### **Features**

- Preview scan (1x to 10x) "Fluorescence" with zoom function
- Preview scan (1x to 10x) "Brightfield" with zoom function
- Automatic multi-channel acquisition of fluorescence images (up to 10 channels)
- Freely selectable magnification for acquisition
- Accurate autofocus for 1x to 100x oil objectives
- One push white balance function
- Detail Windows with zoom in and zoom out function
- Generation of overview images using TissueStitching function
- Export function for detail images and overview with selectable resolution and file size calculation

#### **Benefits**

- Easy to use Project Wizard
- Objective quantification
- Observer independent results
- Easy Handling
- Accuracy
- Freely selectable regions of interest (ROI)
- Selection of individual fields of views by setting flags for re-acquisition or export
- Relocate function
- Customised colours and controlable intensity in overlay images
- Virtual slide creation
- Effective bleaching prevention by minimisation of open shutter times
- Saving of time and resources due to high degree of automatisation
- High throughput acquisition
- Forward and Backward Gating



#### Fast, Precise and Automated Image Acquisition! Easy, observer independent automated tissue analysis!

HAL 100

# **TissueFAXS** i Inverted Cell Analysis System

TissueFAXS i<sup>®</sup> provides the TissueFAXS capabilities on an inverted microscope basis.

This way, quantitative analysis can easily be performed on chamber slides, well plates, microtiter plates, culture flasks as well as standard slides and smears.

The system can be optionally upgraded with third party devices for other inverted microscope applications.

Very detailed information about marker distribution, marker localisation and molecular transport on a sub cellular basis is provided by TissueFAXS i.

Due to this function experiments like intranuclear translocation studies can be easily analysed.

olibri

The analytical capabilities of TissueFAXS i offer the possibility to efficiently implement quality control measures in the pharmaceutical field as well as in cell therapy, where the viability, the phenotype and the activity status of cells can be monitored.



# **Technical Specifications**

**Supported high-end microscopes** 

Zeiss AxioImager.Z1 Zeiss AxioObserver.Z1 Leica DM 6000 Nikon Eclipse 90i

Fully motorized base Offers complete automation. Up to 7 objective lenses From 1x to 100x immersion oil Up to 10 fluorescence filter cubes DAPI, Alexa 488, Cy3, APC are standard

**High-performance workstation** 2 x 24" TFT screens, Intel Core 2 Quad, 4 Gigabyte RAM, 1 Terabyte HD

Required space 160cm x 80cm (5.11ft x 2.55ft)

#### Illumination

Latest fluorescence fiber illumination LED Diascopic illumination (option) Halogen 12V-100W (standard)

#### **High-precision motorized stage**

Stage for upright microscope For up to 8 slides. Stage for invert microscope For microtiter plates and up to 4slides. Resolution Step size as low as 1.5 nm (.0015 μm) for smooth movement. Repeatability Relocation difference < 1 μm, so you can find on the slide precisely what you see on screen. Mark and find Finds absolute stage positions at < 4 μm accu-

racy (as is necessary for reliable automation) Autofocus for brightfield and fluorescence applications

# High-performance monochrome camera for fluorescence

#### Superior quantum efficiency

#### Frame rate

12 fps @ 1392 × 1024 for rapid specimen acquisition.

#### Sensor format

2/3" offering wide fields of view. Resolution

1.4 megapixels guarantees fine object details. Dynamic range

69.5 dB covering different expression levels of biological markers.

#### Acquisition times

Focus	10x len	IS	20x lens		
field	800x	1600x	800x	1600x	
	600	1200	600	1200	
3x3	43s	75s	117s	227s	
5x5	25s	37s	81s	152s	
7x7	26s	37s	67s	123s	

# High-performance colour camera for brightfield (optional)

#### Resolution Up to 1.9 megapixels (1600x1200px) guarantees fine objective detail. Frame rate Up to 100 fps @ VGA for acquisition. Sensor format Up to 1" offers extremely wide fields of view. Excellent colour reproduction Use of multiple chromophores.

# Technical and application support via remote control

#### **TissueGnostics**

# **TissueGnostics** Publication list

Is Benign Prostatic Hyperplasia (BPH) an Immune Inflammatory Disease? Kramer G, Mitteregger D, Marberger M. - Eur Urol. 2007 May 51 (5): 1202-16. Epub 2006 Dec 11.

# Glycogen Synthase Kinase 3beta (GSK3beta) Regulates Differentiation and Proliferation in Neural Stem Cells from the Rat Subventricular Zone Maurer MH, Bromme JO, Feldmann RE Jr, Jarve A, Sabouri F, Burgers HF, Schelshorn DW, Kruger C, Schneider A, Kuschinsky W. J. Proteome - Res. 2007 Mar; 6(3):1198-208.

3D Parallel Coordinate Systems – A New Data Visualization Method in the Context of Microscopy-Based Multicolour Tissue Cytometry (MMTC) M Streit, Rupert C Ecker, K Österreicher, Georg Steiner, H Bischof, C Bangert, T Kopp, R Rogojanu -Cytometry Part

A, 2006 Jul; 69(7):601-11.

An improved method for discrimination of cell populations in tissue sections using Mi-croscopy-Based Multicolour Tissue Cytometry Rupert C Ecker, Radu Rogojanu, Marc Streit, Katja Oesterreicher and Georg Steiner - 2006, Cytometry Part A,

69(3):119-23.

Leukocyte segmentation and classification in blood-smear images Herbert Ramoser, Vincent Laurain, Horst Bischof, Rupert C Ecker - Proc. Of the Int. Conf. Of the IEEE Engineering In Medicine and Biology Society (EMBS); September 01-04, 2005.

Fast Automatic Segmentation of Nuclei in Microscopy Images of Tissue Sections Vincent Laurain, Herbert Ramoser, Christoph Nowak, Rupert C Ecker - Proc. Of the Int. Conf. Of the IEEE Engineer-ing In Medicine and Biology Society (EMBS); September 01-04, 2005.

#### Pimecrolismus leads to an apoptosis-induced depletion of T cells but not Langerhans cells in patients with atopic dermatitis

W Hoetzenecker, Rupert C Ecker, T Kopp, A Stuetz, G Stingl, A Elbe-Bürger - 2005, Journal of Allergy and Clinical Immunology 115(6): 1276-83.

Autocrine IL-10 partially prevents differentiation of neonatal dendritic epidermal leuko-cytes into Langerhans cells

S.Chang-Rodriguez, Rupert C Ecker, G Stingl, A Elbe-Bürger - 2004, Journal of Leukocyte Biology: 76(3): 657-66.

Application of Spectral Imaging Microscopy in Cytomics and Fluorescence Resonance Energy Transfer (FRET) Analysis Rupert C Ecker, R de Martin, GE. Steiner, J A. Schmid - 2004, Cytometry 59A(2): 182-90.

Microscopy-Based Multicolour Tissue Cytometry at the Single Cell Level

Rupert C Ecker and Georg. Steiner - 2004, Cyometry 59A(2):172-81.

Corticosteroids but not Pimecrolimus affect viability, maturation and immune function of murine epidermal Langerhans cells. W Hoetzenecker, JG. Meingassner, Rupert C Ecker, G Stingl, A Stuetz, A Elbe-Buerger - 2004, Journal of Investiga-tive Dermatology 122(3): 673-84.

Inhibition of Restenosis by Tissue Factor Pathway Inhibitor: in vivo and in vitro evidence for suppressed Monocyte Chemoattraction and Gelatinase Activity CW. Kopp, T Hölzenbein, S Steiner, R Marculescu, H Bergmeister, D Seidinger, I Mosberger, C Kaun, M Cejna, R Horvat, J Wojta, G Maurer, Rupert C Ecker, R de Martin, E Minar - 2004, Blood 103(5): 1653-61.

Overexpression of anti-CD75 reactive proteins on distal and collecting renal tubular epi-thelial cells in calcium-oxalate stone forming kidneys in Egypt G Kramer, Georg Steiner, C Neumayer, M Prinz-Kashani, M Hohenfellner, M Gomha, M Ghoneim, M Newman, M Marberger - 2004, British Journal of Urology International 93(6): 822-826.

**Cell-suface matrix proteins and sialic acids in cell-crystal adhesion; the effect of cristal binding on the viability of human CAKI-1 renal epithelial cells** G Kramer, Georg Steiner, M Prinz-Kashani, B Bursa, M Marberger - 2003, British Journal of Urology International 91:

554-559

Published by: TissueGnostics GmbH 1020 Vienna, Austria Tel. +43-1-2161190 Fax +43-1-2161190 90 www.tissuegnostics.com

Automated data acquisition by confocal laser scanning microscopy and image analysis of triple stained immunofluorescent leukocytes in tissue Georg Steiner, Rupert C Ecker, G Kramer, F Stockenhuber, M J Marberger - 2000, J Immunol Methods 237(1-2):



#### Histopathologische Quantifizierung immunologischer Prozesse Rupert C Ecker, R de Martin, and Georg Steiner - 2003, Biomedizinische Technik 48 (Suppl. 1): 546-547.



TissueGnostics GmbH Taborstrasse 10/2/8 A-1020 Wien, Austria Tel.: +43 1 2161190 Fax: +43 1 2161190 90 www.tissuegnostics.com E-mail: office@tissuegnostics.com TissueGnostics USA Ltd. 420 N Larchmont Boulevard CA 90004, Los Angeles, USA Tel.: +1 323 4660499 Fax: +1 323 4660499 www.tissuegnostics.com E-mail: office@tissuegnostics-usa.com TissueGnostics Romania SRL Nicolae lorga Boulevard nr. 51C 700213 Iasi, Romania Tel.: +40 332 405866 Fax: +40 332 405867 www.tissuegnostics.com E-mail: office@tissuegnostics.ro



Copyright<sup>©</sup> TissueGnostics GmbH Taborstrasse10/2/8 - 1020 Vienna - Austria