

# Cover slipper MCS

One of the worldwide fastest automated units for efficient cover slipping



Microtomes

\_\_\_\_\_

Parattin

Embedding Centers

Tissue processors

Stainers

Cover Slipper

► MCS

Knife sharpener

Floating bath

Slide warmer

# Cover slipper MCS

One of the worldwide fastest automated units for efficient cover slipping



The innovative cover slipping technology allows fast identification of single stainings as well as clear classification. The MCS works with highest quality standards, cover slipping is done for each single glass separately fast and without bubbles. This method guarantees a reliable throughput.

The cover slipper MCS can be ideally used in combination with the stainer MAS. The modular design simplifies the step of manual hand over of the sample holders to full automation – also in times of small budgets. The combination uses a minimum of space.

- ► Capacity 120 pieces continuously up to 520 slides per hour
- ► Suitable for commercial cover slipping medium containers
- ► Universal power input 100 240 V, 50/60 Hz
- Active carbon filter for user safety

### **Technical Data**

Cover slipping time	7 – 12 sec. per slide
Cover slipping amount	40 – 250 µl, adjustable

Specimen glasses

Dimensions standard Dimensions optional


Thickness Capacity

## Weight and Dimensions

Dimensions (in mm) Dimensions packed (in mm)

Weight net Weight packed

**Electrical parameters** 



ρ

24 x 40, 50, 55, 60 mm

Depth x Width x Height

24 x 32, 36 mm, 25 x 40, 50, 55, 60 mm

0,12 - 0,17 mm

670 x 900 x 640

800 x 1200 x 800

350 pcs.

75 kg

85 kg

Please ask for other units for Histopathology and consumables as well as laminar flow equipment



### ООО «БиоГен-Аналитика»

115093, Москва, Партийный пер., д.1, корп. 58, стр.1 тел./факс: +7 499 704 62 44 e-mail: 84997046244@bga.su www.bga.su



Stainer MAS in combination with Cover Slipper MCS



Stainer MSM for universal use

### **SLEE MEDICAL GMBH**

Specifications are subject to change without notice)

P. O. Box 10 05 07 D-55129 Mainz · Germany

Fon +49 (0) 6131/95871-0 Fax +49 (0) 6131/95871-0

mail@sleemainz.com www.sleemainz.com