

Rotary Microtome CUT 5062

Precise semi-automation in modular design



Microtomes
CUT 4062
▶ CUT 5062
CUT 6062
MTR

Cryostats

Paraffin

Embedding Centers

Tissue processors

Stainers

Cover slipper

Knife sharpener

Floating bath

Slide warmer

Rotary Microtome CUT 5062

Precise semi-automation in modular design



Standard unit for paraffin sections and research-, plastic- and industrial applications. The semiautomatic unit offers additional advantages against our manual microtomes whilst retaining all essential features of this series. Additional features are:

- ▶ *Motorised trimming up to 300 µm, free programmable*
- ▶ *Motorised advance / return of specimen*
- ▶ *Advance automatic in 2 speeds*
- ▶ *Specimen size up to 50 x 50 mm with standard object clamp or standard cassette size for universal cassette clamp*
- ▶ *Universal power input 100 – 240 V, 50/60 Hz*
- ▶ *Equipped with retraction, free programmable up to 200 µm*
- ▶ *Electric automatic stop function*
- ▶ *The equipment is supplied completely useable including:*
 - *Object orientation unit alternatively fixed*
 - *Universal cassette clamp alternatively standard object clamp*
 - *Knife holder base with disposable blade holder low profile or high profile alternatively standard knife holder*
 - *Section waste tray*



The operating pad for advance / retraction, slow motion, cutting thickness adjustment and trimming



CUT 4062 – manual model



CUT 6062 – fully automatic model

Technical Data

Range of thickness	0,5 µm – 100 µm
Increments	0,5 µm from 0,5 µm – 2 µm 1,0 µm from 2 µm – 10 µm 2,0 µm from 10 µm – 60 µm 5,0 µm from 50 µm – 100 µm
Horizontal specimen advance	28 mm
Vertical stroke	60 mm
Specimen orientation	XY-Axes: 8°, Z-Axes: 360°
Advance for trimming	adjustable from 0 – 300 µm, free programmable

Weight and Dimensions	Depth x Width x Height
Dimensions (in mm)	495 x 400 x 290
Dimensions packed (in mm)	650 x 550 x 480

Weight net	27 kg
Weight packed	39 kg

Electrical parameters	Art.-No. 1070001	230 V / 50 Hz
	Art.-No. 1071001	230 V / 60 Hz
	Art.-No. 1072001	115 V / 60 Hz
	Art.-No. 1073001	100 V / 50 Hz
	Art.-No. 1074001	100 V / 60 Hz

Please ask for other units for Histopathology and consumables as well as laminar flow equipment.

Specifications are subject to change without notice



ООО «БиоГен-Аналитика»

115093, Москва, Партийный пер., д.1, корп. 58, стр.1
 тел./факс: +7 499 704 62 44
 e-mail: 84997046244@bga.su
 www.bga.su

SLEE MEDICAL GMBH

P. O. Box 10 05 07
 D-55129 Mainz · Germany

Fon +49 (0) 6131/95871-0
 Fax +49 (0) 6131/95871-0

mail@sleemainz.com
 www.sleemainz.com