

Tissue Floating Bath MWB

Electrically heated water bath



[Microtomes](#)

[Cryostats](#)

[Paraffin](#)

[Embedding Centers](#)

[Tissue processors](#)

[Stainers](#)

[Cover slipper](#)

[Knife sharpener](#)

[Floating bath](#)

▶ **MWB**

[Slide warmer](#)

Tissue Floating Bath MWB

Electrically heated water bath

SLEE
M A I N Z

- ▶ Bath inside black teflon-coated
- ▶ Heating capacity approx. 65 W
- ▶ Bath with loose stainless steel lid
- ▶ Heating by means of heated foil
- ▶ Double over temperature security integrated in heated circuit
- ▶ Bath fitted in outside wall
- ▶ Outside wall powder coated
- ▶ Electronic temperature regulator with digital temperature display, fitted flush in the front wall
- ▶ CEE-supply socket with main switch and fuses mounted at the rear
- ▶ Temperature selection up to 100°C
- ▶ Comfortable, user friendly spout



Tissue floating bath MWB



Slide warmer MST and floating bath MWB with lid

Technical Data

Bath contents	approx. 3,5 litres
Weight and Dimensions	Depth x Width x Height
Dimensions (in mm)	415 x 275 x 90
Dimensions bath (in mm)	325 x 265 x 65
Dimensions packed (in mm)	500 x 400 x 200
Weight net	3 kg
Weight packed	6 kg
Electrical parameters	Art.-No. 13000200 230 V / 50 Hz Art.-No. 13000201 230 V / 60 Hz Art.-No. 13000202 115 V / 60 Hz Art.-No. 13000203 100 V / 50 Hz Art.-No. 13000204 100 V / 60 Hz

Please ask for other units for Histopathology and consumables as well as laminar flow equipment.

(Specifications are subject to change without notice)



ООО «БиоГен-Аналитика»

115093, Москва, Партийный пер., д.1, корп. 58, стр.1
тел./факс: +7 499 704 62 44
e-mail: 84997046244@bga.su
www.bga.su

SLEE MEDICAL GMBH

P. O. Box 10 05 07
D-55129 Mainz · Germany

Fon +49 (0) 6131/95871-0
Fax +49 (0) 6131/95871-0

mail@sleemainz.com
www.sleemainz.com